

## Antibodies Product Datasheet

|                               |                                                                                                 |
|-------------------------------|-------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | BIOTIN*Monoclonal mouse Anti Human Ig LAMBDA                                                    |
| <b>Catalog Number:</b>        | SA0311                                                                                          |
| <b>Reactivity:</b>            | HumanIgLAMBDA                                                                                   |
| <b>Conjugate:</b>             | Biotin                                                                                          |
| <b>Applications:</b>          | WB,IHC-p,IF/ICC,ELISA                                                                           |
| <b>Dilutions:</b>             | WB 1:5000-20000, ELISA 1:10000-20000, IHC 1:5000-10000 IF 1:5000-10000                          |
| <b>Modification:</b>          | Not Modification Specific                                                                       |
| <b>Source:</b>                | Monoclonal Mouse                                                                                |
| <b>Purification:</b>          | The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen. |
| <b>Concentration:</b>         | 4-5mg/mL                                                                                        |
| <b>Other Name:</b>            | BIOTIN*Monoclonal mouse Anti- Human Ig LAMBDA                                                   |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                         |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                    |