

Antibodies Product Datasheet

Product Name:	BIOTIN*Monoclonal mouse Anti Human IgG4
Catalog Number:	SA0314
Reactivity:	HumanIgG4
Conjugate:	Biotin
Applications:	WB,IHC-p,IF/ICC,ELISA
Dilutions:	WB 1:5000-20000, ELISA 1:10000-20000, IHC 1:5000-10000 IF 1:5000-10000
Modification:	Not Modification Specific
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen.
Concentration:	4-5mg/mL
Other Name:	BIOTIN*Monoclonal mouse Anti- Human IgG4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage and Stability:	-20°C/1 year