

## Antibodies Product Datasheet

|                               |                                                                                                                                                                                                               |
|-------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | Goat Anti Human IgM mu chain HRP                                                                                                                                                                              |
| <b>Catalog Number:</b>        | SA0324                                                                                                                                                                                                        |
| <b>Reactivity:</b>            | Human                                                                                                                                                                                                         |
| <b>Conjugate:</b>             | HRP                                                                                                                                                                                                           |
| <b>Applications:</b>          | ELISA, WB, IHC, ICC                                                                                                                                                                                           |
| <b>Dilutions:</b>             | ELISA (1/5000 - 1/20000), WB (1/5000 - 1/20000), IHC (1/500 - 1/1000), ICC (1/500 - 1/1000)                                                                                                                   |
| <b>Modification:</b>          | Not Modification Specific                                                                                                                                                                                     |
| <b>Source:</b>                | Goat                                                                                                                                                                                                          |
| <b>Purification:</b>          | The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.                                                                                              |
| <b>Concentration:</b>         | 1mg/mL                                                                                                                                                                                                        |
| <b>Other Name:</b>            | Goat Anti-Human IgM mu chain HRP                                                                                                                                                                              |
| <b>Formulation:</b>           | Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide                                                                                                        |
| <b>Storage and Stability:</b> | Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing |