

Antibodies Product Datasheet

Product Name:	Rabbit Anti Guinea Pig IgG (H&L) HRP
Catalog Number:	SA0336
Reactivity:	Guinea Pig
Conjugate:	HRP
Applications:	ELISA, WB, IHC, ICC
Dilutions:	ELISA (1/5000 - 1/20000), WB (1/5000 - 1/20000), IHC (1/500 - 1/1000), ICC (1/500 - 1/1000)
Modification:	Not Modification Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Concentration:	1mg/mL
Other Name:	Rabbit Anti-Guinea Pig IgG (H&L) HRP
Formulation:	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide
Storage and Stability:	Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing