

Antibodies Product Datasheet

Product Name:	DyLight 405, Rabbit Anti Goat IgG
Catalog Number:	SA0394
Reactivity:	Goat
Conjugate:	Dylight 405
Applications:	IF/ICC,FCM
Dilutions:	Optimal working dilutions should be determined experimentally by the investigator. Suggested starting 1:50-1:1000 dilutions for most fluorescent applications.
Modification:	Not Modification Specific
Source:	Polyclonal Rabbit
Purification:	Affinity purified using solid phase Goat IgG (H&L) with finally > 95% purity based on SDS-PAGE
Concentration:	
Other Name:	DyLight 405, Rabbit Anti-Goat IgG
Formulation:	Liquid in PBS, pH 7.4, containing 0.02% Sodium Azide as preservative, 1% BSA as stablizer and 50% Glycerol.
Storage and Stability:	Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing