

Antibodies Product Datasheet

Product Name:	Rabbit Anti Goat IgG (H&L) FITC
Catalog Number:	SA0406
Reactivity:	Goat
Conjugate:	FITC
Applications:	IF/ICC,FCM
Dilutions:	IF (1/250 - 1/2000), FCM (1/250 - 1/2000)
Modification:	Not Modification Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Concentration:	1mg/mL
Other Name:	Rabbit Anti-Goat IgG (H&L) FITC
Formulation:	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide
Storage and Stability:	Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing