

Antibodies Product Datasheet

Product Name:	Mouse Anti-Human Ig kappa chain - AMCA
Catalog Number:	SA0479
Reactivity:	Human
Conjugate:	AMCA
Applications:	ELISA, IHC, IF, FCM
Dilutions:	IHC (1/200 - 1/1000), IF (1/200 - 1/1000), FCM (1/2000 - 1/4000), ELISA (Use at an assay dependent concentration)
Modification:	Not Modification Specific
Source:	Monoclonal Mouse
Purification:	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Concentration:	
Other Name:	
Formulation:	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide
Storage and Stability:	-20°C/1 year