

Antibodies Product Datasheet

Product Name:	Mouse Anti Chicken IgY (H&L) - Alexa Fluor 488
Catalog Number:	SA0685
Reactivity:	Chicken
Conjugate:	Alexa Fluor 488
Applications:	ELISA, IF/ICC, FCM
Dilutions:	IHC (1/200 - 1/1000), IF (1/200 - 1/1000), FCM (1/2000 - 1/4000), ELISA (Use at an assay dependent concentration)
Modification:	Not Modification Specific
Source:	Monoclonal Mouse
Purification:	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Concentration:	1mg/mL
Other Name:	Mouse Anti-Chicken IgY (H&L) - Alexa Fluor 488
Formulation:	1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide
Storage and Stability:	Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing