

Antibodies Product Datasheet

Product Name:	Mouse Anti Human IgM - Alexa Fluor 532
Catalog Number:	SA0733
Reactivity:	Human
Conjugate:	Alexa Fluor 532
Applications:	ELISA,IF/ICC,FCM
Dilutions:	IHC (1/200 - 1/1000), IF (1/200 - 1/1000), FCM (1/2000 - 1/4000), ELISA (Use at an assay dependent concentration)
Modification:	Not Modification Specific
Source:	Monoclonal Mouse
Purification:	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Concentration:	1mg/mL
Other Name:	Mouse Anti-Human IgM - Alexa Fluor 532
Formulation:	1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide
Storage and Stability:	Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing