

## Antibodies Product Datasheet

|                               |                                                                                                                                                                                                               |
|-------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | Goat Anti Human IgG (H&L) - Alexa Fluor 660                                                                                                                                                                   |
| <b>Catalog Number:</b>        | SA0957                                                                                                                                                                                                        |
| <b>Reactivity:</b>            | Human                                                                                                                                                                                                         |
| <b>Conjugate:</b>             | Alexa Fluor 660                                                                                                                                                                                               |
| <b>Applications:</b>          | ELISA,IF/ICC,FCM                                                                                                                                                                                              |
| <b>Dilutions:</b>             | IHC (1/200 - 1/1000), IF (1/200 - 1/1000), FCM (1/2000 - 1/4000), ELISA<br>(Use at an assay dependent concentration)                                                                                          |
| <b>Modification:</b>          | Not Modification Specific                                                                                                                                                                                     |
| <b>Source:</b>                | Goat                                                                                                                                                                                                          |
| <b>Purification:</b>          | The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.                                                                                             |
| <b>Concentration:</b>         | 1mg/mL                                                                                                                                                                                                        |
| <b>Other Name:</b>            | Goat Anti-Human IgG (H&L) - Alexa Fluor 660                                                                                                                                                                   |
| <b>Formulation:</b>           | 1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide                                                                                               |
| <b>Storage and Stability:</b> | Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing |