

Antibodies Product Datasheet

Product Name:	Rabbit Anti Guinea Pig IgG (H&L) - Alexa Fluor 660
Catalog Number:	SA0974
Reactivity:	Guinea Pig
Conjugate:	Alexa Fluor 660
Applications:	ELISA, IF/ICC, FCM
Dilutions:	IHC (1/200 - 1/1000), IF (1/200 - 1/1000), FCM (1/2000 - 1/4000), ELISA (Use at an assay dependent concentration)
Modification:	Not Modification Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Concentration:	1mg/mL
Other Name:	Rabbit Anti-Guinea Pig IgG (H&L) - Alexa Fluor 660
Formulation:	1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide
Storage and Stability:	Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing