

Antibodies Product Datasheet

Product Name:	Goat Anti Human IgG F(ab') ₂ - Alexa Fluor 680
Catalog Number:	SA1040
Reactivity:	Human
Conjugate:	Alexa Fluor 680
Applications:	IF/ICC,FCM,WB
Dilutions:	Suggested starting 1:50-1:1000 dilutions for most fluorescent applications, while 1:5000-20000 dilutions for fluorescent assays with LI-1:COR® Odyssey equipment would be recommended.
Modification:	Not Modification Specific
Source:	Goat
Purification:	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Concentration:	1mg/mL
Other Name:	Goat Anti-Human IgG F(ab') ₂ - Alexa Fluor 680
Formulation:	1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide
Storage and Stability:	Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing