

## Antibodies Product Datasheet

|                               |                                                                                                                      |
|-------------------------------|----------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | Rabbit Anti-Mouse IgG (H&L) - Cy5.5                                                                                  |
| <b>Catalog Number:</b>        | SA1240                                                                                                               |
| <b>Reactivity:</b>            | Mouse                                                                                                                |
| <b>Conjugate:</b>             | Cy5.5                                                                                                                |
| <b>Applications:</b>          | ELISA, IHC, IF, FCM                                                                                                  |
| <b>Dilutions:</b>             | IHC (1/200 - 1/1000), IF (1/200 - 1/1000), FCM (1/2000 - 1/4000), ELISA<br>(Use at an assay dependent concentration) |
| <b>Modification:</b>          | Not Modification Specific                                                                                            |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                    |
| <b>Purification:</b>          | The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.    |
| <b>Concentration:</b>         |                                                                                                                      |
| <b>Other Name:</b>            |                                                                                                                      |
| <b>Formulation:</b>           | Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide               |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                         |