

Antibodies Product Datasheet

Product Name:	Rabbit Anti-Human IgG F(ab') ₂ - Cy5.5
Catalog Number:	SA1256
Reactivity:	Human
Conjugate:	Cy5.5
Applications:	ELISA, IHC, IF, FCM
Dilutions:	IHC (1/200 - 1/1000), IF (1/200 - 1/1000), FCM (1/2000 - 1/4000), ELISA (Use at an assay dependent concentration)
Modification:	Not Modification Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Concentration:	
Other Name:	
Formulation:	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide
Storage and Stability:	-20°C/1 year