

Antibodies Product Datasheet

Product Name:	Mouse Anti Monkey IgG (H&L) HRP
Catalog Number:	SA1301
Reactivity:	Monkey
Applications:	Elisa, WB, IHC, ICC
Dilutions:	Elisa (1:5000 - 1:20000), WB (1:5000 - 1:20000), IHC (1:500 - 1:1000), ICC (1:500 - 1:1000)
Modification:	Unmodified
Source:	Monoclonal, Mouse
Purification:	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Immunogen:	
Concentration:	1mg/mL
Gene Name:	
Protein Name:	Monkey IgG (H&L) HRP
Other Name:	Mouse Anti-Monkey IgG (H&L) HRP
Human Swiss Prot No.:	
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	
Storage and Stability:	Store at -20° C/1 year