

CD38 (AB1408) mouse mAb

V0496

Key Features

Host Species

- Mouse

Reactivity

- Human;

Applications

- Immunohistochemistry (IHC)

MW

- 34 kDa (calculated)
- 40 kDa (observed)

Isotype

- IgG2b, kappa

Recommended Dilution Ratios

Application

Immunohistochemistry (IHC)

Dilution

IHC, 1:200-1:400

Storage

Storage Conditions

Store at -20°C. Avoid freeze / thaw cycles.

Storage buffer

The antibody is provided in liquid form in phosphate - buffered saline with 50% glycerol, 0.05% BSA, and 0.05% Proclin 300.

Basic Information

Clonality Monoclonal

Clone Number AB1408

Immunogen A synthetic peptide corresponding to the amino acid region 200 - 300 of the human CD38 protein.

Specificity The antibody can specifically recognize human CD38 protein.

Purification Affinity purification Protein A

Concentration Product concentration may vary by batch. Please refer to the product COA for details.

Target Information

Gene name CD38

Protein Name CD38

Database Link	Organism	Swiss Prot.	Gene ID
	Human	P28907	952

Background CD38 is a type II transmembrane glycoprotein with a molecular weight of 42 kDa. It is mainly expressed in plasma cells, pre-B cells, activated B cells, thymocytes and activated T cells. 90% of CD34 positive cells also express CD38. It is used to study the role of activated T cells in autoimmune diseases and immune deficiency diseases, and to diagnose multiple myeloma (plasmacytoma) . It is a useful marker in the immunotyping of acute leukemia.