

Antibodies Product Datasheet

Product Name:	Nrf2 (PTR2557) mouse mAb
Catalog Number:	W0037
Reactivity:	Human, Mouse, Rat
Applications:	WB; IF; ELISA
Dilutions:	WB 1:500-2000. IF 1:100-500. ELISA 1:1000-5000
Modification:	Unmodified/Total
Source:	Mouse, Monoclonal/IgG
Purification:	Protein G
Immunogen:	Synthesized peptide derived from human Nrf2 AA range: 300-400
Concentration:	
Gene Name:	NFE2L2 NRF2
Protein Name:	Nuclear factor erythroid 2-related factor 2 (NF-E2-related factor 2) (NFE2-related factor 2) (HEBP1) (Nuclear factor, erythroid derived 2, like 2)
Other Name:	Nuclear factor erythroid 2-related factor 2 (NF-E2-related factor 2) (NFE2-related factor 2) (HEBP1) (Nuclear factor, erythroid derived 2, like 2)
Human Swiss Prot No.:	Q16236
Formulation:	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Specificity:	This antibody detects endogenous levels of Nrf2 protein.
Storage and Stability:	Store at -20° C/1 year