

Antibodies Product Datasheet

Product Name:	eIF2 α (Phospho Ser51) (PTR1316) mouse mAb
Catalog Number:	W0159
Reactivity:	Human, Mouse, Rat
Applications:	WB, IF, ELISA
Dilutions:	
Modification:	Unmodified/Total
Source:	
Purification:	
Immunogen:	
Concentration:	
Gene Name:	EIF2S1
Protein Name:	Eukaryotic translation initiation factor 2 subunit 1; Eukaryotic translation initiation factor 2 subunit alpha (eIF-2-alpha; eIF-2A; eIF-2 α); EIF2S1; EIF2A
Other Name:	
Human Swiss Prot No.:	P05198
Formulation:	
Specificity:	
Storage and Stability:	Store at -20° C/1 year