

Antibodies Product Datasheet

Product Name:	Glial Fibrillary Acidic Protein (GFAP) (PTR2201/470) mouse mAb
Catalog Number:	W0188
Reactivity:	Human, Mouse, Rat, Monkey, Bovine,
Applications:	IHC; IF; ELISA
Dilutions:	IHC 1:200-1000. IF 1:100-500. ELISA 1:1000-5000
Modification:	Unmodified/Total
Source:	Mouse, Monoclonal/IgG1, kappa
Purification:	Protein G
Immunogen:	Synthesized peptide derived from human GFAP AA range: 300-432
Concentration:	
Gene Name:	GFAP
Protein Name:	GFAP
Other Name:	wu:fb34h11 ALXDRD cb345 etID36982.3 FLJ42474 FLJ45472 GFAP GFAP_HUMAN gfapl Glial fibrillary acidic protein Intermediate filament protein wu:fk42c12 xx:af506734 zgc:110485
Human Swiss Prot No.:	P14136
Formulation:	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Specificity:	This antibody detects endogenous levels of GFAP protein.
Storage and Stability:	Store at -20° C/1 year