

Antibodies Product Datasheet

Product Name:	Lamin A/C (phospho Ser392) (PTR1133) mouse mAb
Catalog Number:	W0215
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, IF, ELISA
Dilutions:	
Modification:	Unmodified/Total
Source:	
Purification:	
Immunogen:	
Concentration:	
Gene Name:	LMNA LMN1
Protein Name:	Prelamin-A/C [Cleaved into: Lamin-A/C (70 kDa lamin) (Renal carcinoma antigen NY-REN-32)]
Other Name:	
Human Swiss Prot No.:	P02545
Formulation:	
Specificity:	
Storage and Stability:	Store at -20° C/1 year