

CD8 a (AB2392) mouse mAb

W0232

Key Features

Host Species

- Mouse

Reactivity

- Mouse;

Applications

- IHC; WB; IF; ELISA

MW

- 26 kDa (calculated)
- 30-35 kDa (observed)

Isotype

- IgG2b, kappa

Recommended Dilution Ratios

Application

IHC; WB; IF; ELISA

Dilution

IHC, 1:200-1:1000 | WB, 1:500-1:2000 | IF, 1:100-1:500 | ELISA, Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Storage

Storage Conditions

Store at -20°C. Avoid freeze / thaw cycles.

Storage buffer

The antibody is provided in liquid form in phosphate - buffered saline with 50% glycerol, 0.05% BSA, and 0.05% Proclin 300.

Basic Information

Clonality	Monoclonal
Clone Number	AB2392
Immunogen	A synthetic peptide corresponding to the amino acid region 150 - 235 of the mouse CD8 a protein.
Specificity	This antibody detects endogenous levels of mouse CD8 a protein.
Purification	Affinity purification Protein A
Concentration	Product concentration may vary by batch. Please refer to the product COA for details.

Target Information

Gene name	CD8A MAL		
Protein Name	CD8 a		
Database Link	Organism	Swiss Prot.	Gene ID
	Human	P01732	925

Background

The CD8 antigen is a cell surface glycoprotein found on most cytotoxic T lymphocytes that mediates efficient cell-cell interactions within the immune system. The CD8 antigen acts as a coreceptor with the T-cell receptor on the T lymphocyte to recognize antigens displayed by an antigen presenting cell in the context of class I MHC molecules. The coreceptor functions as either a homodimer composed of two alpha chains or as a heterodimer composed of one alpha and one beta chain. Both alpha and beta chains share significant homology to immunoglobulin variable light chains. This gene encodes the CD8 alpha chain. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2011],



2200 Ringwood Ave. San Jose, CA 95131
1 (877) 594-3616 (Toll Free), (408) 747-0185

www.assaybiotechnology.com
tech@assaybiotech.com | order@assaybiotech.com
