

Claudin 3 (AB1534) mouse mAb

W0329

Key Features

Host Species

- Mouse

Reactivity

- Human;

Applications

- IHC; WB; IF; ELISA

MW

- 23 kDa (calculated)
- 20 kDa (observed)

Isotype

- IgG1, kappa

Recommended Dilution Ratios

Application

IHC; WB; IF; ELISA

Dilution

IHC, 1:200-1:1000 | WB, 1:500-1:2000 | IF, 1:100-1:500 | ELISA, Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Storage

Storage Conditions

Store at -20°C. Avoid freeze / thaw cycles.

Storage buffer

The antibody is provided in liquid form in phosphate - buffered saline with 50% glycerol, 0.05% BSA, and 0.05% Proclin 300.

Basic Information

Clonality	Monoclonal
Clone Number	AB1534
Immunogen	A synthetic peptide corresponding to the amino acid region 150 - 220 of the human Claudin 3 protein.
Specificity	This antibody detects endogenous levels of Claudin 3 protein.
Purification	Affinity purification Protein A
Concentration	Product concentration may vary by batch. Please refer to the product COA for details.

Target Information

Gene name	CLDN3		
Protein Name	Claudin 3		
Database Link	Organism	Swiss Prot.	Gene ID
	Human	O15551	1365

Background

Tight junctions represent one mode of cell-to-cell adhesion in epithelial or endothelial cell sheets, forming continuous seals around cells and serving as a physical barrier to prevent solutes and water from passing freely through the paracellular space. These junctions are comprised of sets of continuous networking strands in the outwardly facing cytoplasmic leaflet, with complementary grooves in the inwardly facing extracytoplasmic leaflet. The protein encoded by this intronless gene, a member of the claudin family, is an integral membrane protein and a component of tight junction strands. It is also a low-affinity receptor for *Clostridium perfringens* enterotoxin, and shares aa sequence similarity with a putative apoptosis-related protein found in rat. [provided by RefSeq, Jul 2008]



2200 Ringwood Ave. San Jose, CA 95131
1 (877) 594-3616 (Toll Free), (408) 747-0185

www.assaybiotechnology.com
tech@assaybiotech.com | order@assaybiotech.com
