

E-Cadherin (AB1337) mouse mAb

W0340

Key Features

Host Species

- Mouse

Reactivity

- Human; Mouse; Rat;

Applications

- IHC; WB; IF; ELISA

MW

- 97 kDa (calculated)
- 120 kDa (observed)

Isotype

- IgG2b, kappa

Recommended Dilution Ratios

Application

IHC; WB; IF; ELISA

Dilution

IHC, 1:200-1:1000 | WB, 1:500-1:2000 | IF, 1:100-1:500 | ELISA, Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Storage

Storage Conditions

Store at -20°C. Avoid freeze / thaw cycles.

Storage buffer

The antibody is provided in liquid form in phosphate - buffered saline with 50% glycerol, 0.05% BSA, and 0.05% Proclin 300.

Basic Information

Clonality	Monoclonal
Clone Number	AB1337
Immunogen	A synthetic peptide corresponding to the amino acid region 700 - 800 of the human E - Cadherin protein.
Specificity	The antibody can specifically recognize E-cadherin protein.
Purification	Affinity purification Protein A
Concentration	Product concentration may vary by batch. Please refer to the product COA for details.

Target Information

Gene name	CDH1		
Protein Name	E-Cadherin		
Database Link	Organism	Swiss Prot.	Gene ID
	Human	P12830	999

Background

E-Cadherin, 120kDa, is a critical regulator of epithelial junction formation. The E-cadherin internal domain binds with alpha, beta, gamma, and p120 catenins to anchor the ECAD complex to the actin cytoskeleton of the cell. E-cadherin is expressed in virtually all epithelial cells with the exception of adrenocortical cells. The expression in liver cells is weaker than in most other epithelia. ECAD is also expressed in melanocytes (adhering to squamous epithelial cells) . Mesothelial cells and granulosa cells are ECAD negative. Mesenchymal cells and neural crest cells are mostly negative but some types may show a weak reaction. Testicular germinative cells are weakly positive. At present, it is commonly used to differentiate breast lobular cancer from ductal cancer (the former is often negative, the latter is often positive) .



2200 Ringwood Ave. San Jose, CA 95131
1 (877) 594-3616 (Toll Free), (408) 747-0185

www.assaybiotechnology.com
tech@assaybiotech.com | order@assaybiotech.com
