

## Assay Biotechnology Company, Inc.

1276 Hammerwood Ave.

Sunnyvale, CA 94089

Phone: (408) 747-0185

Fax: (408)747-0145

www.assaybiotech.com

## NF-kB p65 (Phospho-Ser536) Phospho Sandwich ELISA Kit

Catalog Number: PS-1633

### Description

#### Details

Product Name: NF-kB p65 (Phospho-Ser536) Phospho Sandwich ELISA Kit

Catalog No: PS-1633

Stock:

#### Properties

Dynamic Range: >16 µg/mL Total Protein

Storage/Stability: 4°C/6 Months

Detection Method: Colorimetric 450 nm

Reactivity: Human, Mouse, Rat

Assay Type: Sandwich ELISA

Database Links: Gene ID: 5970, UniProt ID: Q04206

Format: 12 x 8-Well Microstrips

Image Legend: HeLa cells were grown to 90% confluency and were stimulated with EGF 200ng/mL for 30 minutes. Cells were immediately lysed thereafter and measured for Total Protein Concentration and O.D. 450 nm of NFkB P-536 versus untreated HeLa lysates.

Manual (PDF): [Download](#)

NCBI Gene  
Symbol: RELA

Specificity: NF-kB p65 (Phospho-Ser536) Sandwich ELISA Kit detects endogenous levels of NF-kB p65 (Phospho-Ser536) protein.

Sub Type: Phosphorylation

Synonyms: NFKB3; nuclear factor NF-kappa-B p65 subunit; p65; RELA; TF65; transcription factor p65

TargetName: NF-kB p65 (Phospho-Ser536)