

## Assay Biotechnology Company, Inc.

1276 Hammerwood Ave.

Sunnyvale, CA 94089

Phone: (408) 747-0185

Fax: (408)747-0145

www.assaybiotech.com

## GSK3 beta (Phospho-Ser9) Phospho Sandwich ELISA Kit

Catalog Number: PS-1567

### Description

#### Details

Product Name: GSK3 beta (Phospho-Ser9) Phospho Sandwich ELISA Kit

Catalog No: PS-1567

Stock:

#### Properties

Dynamic Range: >16 µg/mL Total Protein

Storage/Stability: 4°C/6 Months

Detection Method: Colorimetric 450 nm

Reactivity: Human, Mouse, Rat

Assay Type: Sandwich ELISA

Database Links: Gene ID: 2932, UniProt ID: P49841

Format: 12 x 8-Well Microstrips

Image Legend: HeLa cells were grown to 90% confluency and were stimulated with EGF 200 ng/mL for 30 minutes. Cells were immediately lysed thereafter and measured for Total Protein Concentration and O.D. 450 nm of GSK3 beta P-Ser9 versus untreated HeLa lysates.

Manual (PDF): [Download](#)

NCBI Gene Symbol: GSK3B

Specificity: GSK3 beta (Phospho-Ser9) Sandwich ELISA Kit detects endogenous levels of GSK3 beta (Phospho-Ser9) protein.

Sub Type: Phosphorylation

**Synonyms:** Factor A; GSK-3 alpha; GSK3A; Glycogen synthase kinase-3 alpha; kinase  
GSK3-alpha

**TargetName:** GSK3 beta (Phospho-Ser9)