

Assay Biotechnology Company, Inc.

1276 Hammerwood Ave.

Sunnyvale, CA 94089

Phone: (408) 747-0185

Fax: (408)747-0145

www.assaybiotech.com

Anti-S100A4 antibody [ABT-S100A4]

Catalog Number: V0122

Description

Details

Product Name: Anti-S100A4 antibody [ABT-S100A4]

Catalog No: V0122

Stock:

Properties

Description: Mouse monoclonal [ABT-S100A4] antibody to S100A4

Applications: IHC

Clonality: Monoclonal

Host Species: Mouse

Reactivity: Human

ABT Clone ID: ABT-S100A4

Antigen Retrieval: Heat-induced epitope retrieval (HIER)

Concentration: 1mg/ml

Dilution: 1:50-200

Entrez Gene ID
(NCBI): 6275

IHC
Recommended
control: Malignant melanoma

IHC Target Name: S100A4

Immunogen:	Synthetic peptide
Localization:	Nuclear, Cytoplasmic
OMIM:	114210
Purification:	Immunogen affinity purified
Retrieval buffer:	Citrate buffer of pH6.0
Source/Ig Isotype:	Mouse IgG
Storage/Stablility:	-20°C/1 year
Synonyms:	Protein S100-A4, calvasculin, metastasin, placental calcium-binding protein, protein Mts1, S100 calcium-binding protein A4, CAPL, MTS1
UniProt Gene Name:	S100A4
UniProtKB:	P26447
Volume:	2ml