

Assay Biotechnology Company, Inc.

1276 Hammerwood Ave.

Sunnyvale, CA 94089

Phone: (408) 747-0185

Fax: (408)747-0145

www.assaybiotech.com

Anti-TGF β Receptor II antibody [ABT-TGFR2]

Catalog Number: V0130

Description

Details

Product Name: Anti-TGF β Receptor II antibody [ABT-TGFR2]

Catalog No: V0130

Stock:

Properties

Description: Mouse monoclonal [ABT-TGFR2] antibody to TGFR2

Applications: IHC

Clonality: Monoclonal

Host Species: Mouse

Reactivity: Human

ABT Clone ID: ABT-TGFR2

Antigen Retrieval: Heat-induced epitope retrieval (HIER)

Dilution: 1:50-200

Entrez Gene ID
(NCBI): 7048

IHC
Recommended
control: Colon

IHC Target Name: TGF β Receptor II

Immunogen: Synthetic peptide

| | |
|---------------------|---|
| Localization: | Membranous |
| OMIM: | 190182 |
| Purification: | Immunogen affinity purified |
| Retrieval buffer: | Citrate buffer of pH6.0 |
| Source/Ig Isotype: | Mouse IgG |
| Storage/Stablility: | -20°C/1 year |
| Synonyms: | TGF-beta receptor type II, TGFR-2, TGF-beta type II receptor, Transforming growth factor-beta receptor type II, TbetaR-II |
| UniProt Gene Name: | TGFBR2 |
| UniProtKB: | P37173 |