

Assay Biotechnology Company, Inc.

1276 Hammerwood Ave.

Sunnyvale, CA 94089

Phone: (408) 747-0185

Fax: (408)747-0145

www.assaybiotech.com

FKHRL1 (Phospho-Ser253) Fluorometric Cell-Based ELISA Kit

Catalog Number: FLUO-CBP1558

Description

Details

Product Name: FKHRL1 (Phospho-Ser253) Fluorometric Cell-Based ELISA Kit

Catalog No: FLUO-CBP1558

Stock:

Properties

Detection Method: Fluorometric

Dynamic Range: > 5000 Cells

Storage/Stability: 4°C/6 Months

Reactivity: Human, Mouse, Rat

Assay Type: Cell-Based ELISA

Database Links: Gene ID: 2309, UniProt ID: O43524, OMIM #: 602681, Unigene #: Hs.220950

Format: Two 96-Well Plates

Manual (PDF): [Download](#)

NCBI Gene Symbol: FOXO3

Sub Type: Phospho

Synonyms: AF6q21 protein, FKHR2, FKHRL1, Forkhead box O3, forkhead box O3A, Forkhead box protein O3A, Forkhead in rhabdomyosarcoma-like 1, FOXO3

TargetName: Phospho-FKHRL1 (Ser253)