

Assay Biotechnology Company, Inc.

1276 Hammerwood Ave.

Sunnyvale, CA 94089

Phone: (408) 747-0185

Fax: (408)747-0145

www.assaybiotech.com

53BP1 (Phospho-Ser6) Colorimetric Cell-Based ELISA Kit

Catalog Number: CBP1286

Description

Details

Product Name: 53BP1 (Phospho-Ser6) Colorimetric Cell-Based ELISA Kit

Catalog No: CBP1286

Stock:

Properties

Detection Method: Colorimetric 450 nm

Dynamic Range: > 5000 Cells

Storage/Stability: 4°C/6 Months

Reactivity: Human, Mouse, Rat

Assay Type: Cell-Based

Database Links: Gene ID: 7158, UniProt ID: Q12888, OMIM #: 605230, Unigene #: Hs.440968

Format: 2 x 96-Well Microplates

Manual (PDF): [Download](#)

NCBI Gene Symbol: TP53BP1

Sub Type: Phospho

Synonyms: p53-binding protein 1, P531, TP53B, TP53BP1, TRP53BP1, Tumor suppressor p53-binding protein 1

TargetName: 53BP1 (Phospho-Ser6)