

Assay Biotechnology Company, Inc.

1276 Hammerwood Ave.

Sunnyvale, CA 94089

Phone: (408) 747-0185

Fax: (408)747-0145

www.assaybiotech.com

Phospho-GAB1 (Tyr659) Colorimetric Cell-Based ELISA Kit

Catalog Number: CBP1301

Description

Details

Product Name: Phospho-GAB1 (Tyr659) Colorimetric Cell-Based ELISA Kit

Catalog No: CBP1301

Stock:

Properties

Detection Method: Colorimetric 450 nm

Dynamic Range: > 5000 Cells

Storage/Stability: 4°C/6 Months

Reactivity: Human, Mouse, Rat

Assay Type: Cell-Based

Database Links: Gene ID: 2549, UniProt ID: Q13480, OMIM #: 604439, Unigene #: Hs.618456/Hs.80720

Format: 2 x 96-Well Microplates

Manual (PDF): [Download](#)

NCBI Gene Symbol: GAB1

Sub Type: Phospho

Synonyms: GRB2-associated binder 1, GRB2-associated-binding protein 1, growth factor receptor bound protein 2-associated protein 1

TargetName: Phospho-GAB1 (Tyr659)