

Assay Biotechnology Company, Inc.

1276 Hammerwood Ave.

Sunnyvale, CA 94089

Phone: (408) 747-0185

Fax: (408)747-0145

www.assaybiotech.com

EGF (Human) LumiAb[™]; ELISA Kit

Catalog Number: Lum-8111

Description

Details

Product Name: EGF (Human) LumiAb[™]; ELISA Kit

Catalog No: Lum-8111

Stock:

Properties

Detection Method: Luminometer

Dynamic Range: 4-250

Storage/Stability: 4°C/6 Months

Reactivity: Human

Assay Type: Sandwich ELISA

Database Links: Gene ID: 1950, UniProt ID: P01133

Format: 12 x 8-Well Microstrips

Image Legend: The LumiAb[™], Human EGF ELISA Kit allows for the detection and quantification of endogenous levels of natural and/or recombinant Human EGF proteins within the range of 4-250 pg/ml.

Manual (PDF): [Download](#)

NCBI Gene Symbol: EGF

Specificity: The LumiAb[™], Human EGF ELISA is capable of recognizing both recombinant and naturally produced Human EGF proteins. The antigens listed below were tested at 50 ng/ml and did not exhibit significant cross reactivity or interference. □ Human: EGF-Receptor, FGF-Acidic, FGF-Basic, FGF-9, FGF-18, IL-17, TNF-RI, 4-1BBL, sCD40-L □ Murine: EGF, GM-CSF, TNF-alpha □ Rat: EGF,

TNF-alpha

Sub Type: None

Synonyms: Epidermal growth factor, Urogastrone

TargetName: Human EGF