

Assay Biotechnology Company, Inc.

1276 Hammerwood Ave.

Sunnyvale, CA 94089

Phone: (408) 747-0185

Fax: (408)747-0145

www.assaybiotech.com

IRF3 Antibody

Catalog Number: B0496

Description

Details

Product Name: IRF3 Antibody

Catalog No: B0496

Stock:

Properties

Storage/Stability: Stable at -20°C for at least 1 year.

Applications: WB, IHC, IF, ELISA

Clonality: Polyclonal

Host Species: Rabbit

Reactivity: Human

Application Dilutions: WB: 1:500~1:1000 IHC: 1:50~1:100 IF: 1:100~1:500 ELISA: 1:5000

Database Links: Gene ID: 3661, Swiss-Prot #: Q14653, OMIM #: 603734, Unigene #: Hs.75254

Immunogen: The antiserum was produced against synthesized peptide derived from human IRF3.

Immunogen Range: 351-400

MW: 47 kDa

NCBI Gene Symbol: IRF3

Physical Form: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Purification: The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.

Specificity: IRF3 Antibody detects endogenous levels of total IRF3 protein.

Synonyms: Interferon regulatory factor 3, IRF3

Target Modification: Unmodified/Total