

## Assay Biotechnology Company, Inc.

1276 Hammerwood Ave.  
Sunnyvale, CA 94089  
Phone: (408) 747-0185  
Fax: (408)747-0145  
www.assaybiotech.com

## IRF3 Antibody

Catalog Number: B0667

## Description

### Details

Product Name: IRF3 Antibody

Catalog No: B0667

Stock:

### Properties

Storage/Stability: Stable at -20°C for at least 1 year.

Applications: WB, IHC, ELISA

Clonality: Polyclonal

Host Species: Rabbit

Reactivity: Human, Mouse, Rat

Application Dilutions: WB: 1:500~1:1000 IHC: 1:50~1:100 ELISA: 1:1000

Database Links: Gene ID: 3661, Swiss-Prot #: Q14653, OMIM #: 603734, Unigene #: Hs.75254

Immunogen: The antiserum was produced against synthesized peptide derived from human IRF3.

Immunogen Range: 362-411

MW: 47 kDa

NCBI Gene Symbol: IRF3

Physical Form: Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Purification:** The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.

**Specificity:** IRF3 Antibody detects endogenous levels of total IRF3 protein.

**Synonyms:** Interferon regulatory factor 3, IRF3

**Target Modification:** Unmodified/Total