

## Assay Biotechnology Company, Inc.

1276 Hammerwood Ave.

Sunnyvale, CA 94089

Phone: (408) 747-0185

Fax: (408)747-0145

www.assaybiotech.com

## Anti-Actin pan antibody [ABT-ACTN]

Catalog Number: V0005

### Description

#### Details

Product Name: Anti-Actin pan antibody [ABT-ACTN]

Catalog No: V0005

Stock:

#### Properties

Description: Mouse monoclonal [ABT-ACTN] antibody to ACTN

Applications: IHC, WB

Clonality: Monoclonal

Host Species: Mouse

Reactivity: Human

ABT Clone ID: ABT-ACTN

Antigen Retrieval: Heat-induced epitope retrieval (HIER)

Concentration: 0.6mg/ml

Dilution: 1:50-200

IHC Recommended control: Normal Colon, Cardiac muscle, Skeletal muscle

IHC Target Name: Actin pan

Immunogen: Synthetic peptide

Localization: Cytoplasmic

|                     |                             |
|---------------------|-----------------------------|
| Purification:       | Immunogen affinity purified |
| Retrieval buffer:   | Citrate buffer of pH6.0     |
| Source/Ig Isotype:  | Mouse IgG                   |
| Storage/Stablility: | -20°C/1 year                |
| UniProtKB:          | PAN                         |
| Volume:             | 2ml                         |