

## Assay Biotechnology Company, Inc.

1276 Hammerwood Ave.

Sunnyvale, CA 94089

Phone: (408) 747-0185

Fax: (408)747-0145

www.assaybiotech.com

## Gab2 Antibody

Catalog Number: B8349

## Description

### Details

Product Name: Gab2 Antibody

Catalog No: B8349

Stock:

### Properties

Storage/Stability: Stable at -20°C for at least 1 year.

Applications: WB, ELISA

Clonality: Polyclonal

Host Species: Rabbit

Reactivity: Human, Mouse, Rat

Application Dilutions: WB: 1:500~1:1000 ELISA: 1:40000

Database Links: Gene ID: 9846, Swiss-Prot #: Q9UQC2, OMIM #: 606203, Unigene #: Hs.429434

Immunogen: The antiserum was produced against synthesized peptide derived from human Gab2.

Immunogen Range: 609-658

MW: 74 kDa

NCBI Gene Symbol: GAB2

Physical Form: Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Purification:** The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.

**Specificity:** Gab2 Antibody detects endogenous levels of total Gab2 protein.

**Synonyms:** GRB2-associated-binding protein 2, Growth factor receptor bound protein 2-associated protein 2, GRB2-associated binder 2, pp100, KIAA0571

**Target Modification:** Unmodified/Total