

Assay Biotechnology Company, Inc.

1276 Hammerwood Ave.

Sunnyvale, CA 94089

Phone: (408) 747-0185

Fax: (408)747-0145

www.assaybiotech.com

GRP94 Antibody

Catalog Number: C0218

Description

Details

Product Name: GRP94 Antibody

Catalog No: C0218

Stock:

Properties

Storage/Stability: Stable at -20°C for at least 1 year.

Applications: WB, IHC, IF, ELISA

Clonality: Polyclonal

Host Species: Rabbit

Reactivity: Human, Mouse

Application Dilutions: WB: 1:500~1:1000 IHC: 1:50~1:100 IF: 1:100~1:500 ELISA: 1:5000

Database Links: Gene ID: 7184, Swiss-Prot #: P14625, OMIM #: 191175, Unigene #: Hs.192374

Immunogen: The antiserum was produced against synthesized peptide derived from human GRP94.

Immunogen Range: 754-803

MW: 92 kDa

NCBI Gene Symbol: HSP90B1

Physical Form: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Purification: The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.

Specificity: GRP94 Antibody detects endogenous levels of total GRP94 protein.

Synonyms: TRA1, GLUCOSE-REGULATED PROTEIN, 94-KD, STRESS-INDUCIBLE TUMOR REJECTION ANTIGEN GP96, TUMOR REJECTION ANTIGEN 1, tumor rejection antigen (gp96) 1

Target Modification: Unmodified/Total