

Assay Biotechnology Company, Inc.

1276 Hammerwood Ave.

Sunnyvale, CA 94089

Phone: (408) 747-0185

Fax: (408)747-0145

www.assaybiotech.com

Anti-cadherin pan antibody [ABT-MCDH]

Catalog Number: V0011

Description

Details

Product Name: Anti-cadherin pan antibody [ABT-MCDH]

Catalog No: V0011

Stock:

Properties

Description: Mouse monoclonal [ABT-MCDH] antibody to MCDH

Applications: IHC, WB

Clonality: Monoclonal

Host Species: Mouse

Reactivity: Human

ABT Clone ID: ABT-MCDH

Antigen Retrieval: Heat-induced epitope retrieval (HIER)

Dilution: 1:50-200

IHC
Recommended
control: Colon

IHC Target Name: cadherin pan

Immunogen: Synthetic peptide

Localization: Membranous

Purification: Immunogen affinity purified

Retrieval buffer:	Citrate buffer of pH6.0
Source/Ig Isotype:	Mouse IgG
Storage/Stablility:	-20°C/1 year
UniProtKB:	MANY