

Assay Biotechnology Company, Inc.

1276 Hammerwood Ave.

Sunnyvale, CA 94089

Phone: (408) 747-0185

Fax: (408)747-0145

www.assaybiotech.com

KAPCB Antibody

Catalog Number: C10453

Description

Details

Product Name: KAPCB Antibody

Catalog No: C10453

Stock:

Properties

Storage/Stability: Stable at -20°C for at least 1 year.

Applications: WB, ELISA

Clonality: Polyclonal

Host Species: Rabbit

Reactivity: Human, Mouse, Rat

Application Dilutions: WB: 1:500~1:1000 ELISA: 1:20000

Database Links: Gene ID: 5567, Swiss-Prot #: P22694, OMIM #: 176892, Unigene #: Hs.487325

Immunogen: The antiserum was produced against synthesized peptide derived from human KAPCB.

Immunogen Range: 291-340

MW: 40 kDa

NCBI Gene Symbol: PRKACB

Physical Form: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Purification: The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.

Specificity: KAPCB Antibody detects endogenous levels of total KAPCB protein.

Synonyms: KAPB, KAPCB, PKA C-beta, PKA-beta, cAMP-dependent protein kinase, beta-catalytic subunit

Target Modification: Unmodified/Total