

Assay Biotechnology Company, Inc.

1276 Hammerwood Ave.

Sunnyvale, CA 94089

Phone: (408) 747-0185

Fax: (408)747-0145

www.assaybiotech.com

TCF7 Antibody

Catalog Number: C10571

Description

Details

Product Name: TCF7 Antibody

Catalog No: C10571

Stock:

Properties

Storage/Stability: Stable at -20°C for at least 1 year.

Applications: WB, ELISA

Clonality: Polyclonal

Host Species: Rabbit

Reactivity: Human, Mouse

Application Dilutions: WB: 1:500~1:1000 ELISA: 1:10000

Database Links: Gene ID: 6932, Swiss-Prot #: P36402, OMIM #: 189908, Unigene #: Hs.573153

Immunogen: The antiserum was produced against synthesized peptide derived from human TCF7.

Immunogen Range: 8-59

MW: 41 kDa

NCBI Gene Symbol: TCF7

Physical Form: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Purification: The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.

Specificity: TCF7 Antibody detects endogenous levels of total TCF7 protein.

Synonyms: Transcription factor 7, T-cell-specific transcription factor 1, T-cell factor 1, TCF-1, TCF7, TCF1

Target Modification: Unmodified/Total