

## Assay Biotechnology Company, Inc.

1276 Hammerwood Ave.

Sunnyvale, CA 94089

Phone: (408) 747-0185

Fax: (408)747-0145

www.assaybiotech.com

## Eotaxin-2 (Human) LumiAb<sup>™</sup>; ELISA Kit

Catalog Number: Lum-8309

### Description

#### Details

Product Name: Eotaxin-2 (Human) LumiAb<sup>™</sup>; ELISA Kit

Catalog No: Lum-8309

Stock:

#### Properties

Detection Method: Luminometer

Dynamic Range: 8-500

Storage/Stability: 4°C/6 Months

Reactivity: Human

Assay Type: Sandwich ELISA

Database Links: Gene ID: 6369, UniProt ID: O00175

Format: 12 x 8-Well Microstrips

Image Legend: The LumiAb<sup>™</sup>, Human Eotaxin-2 ELISA Kit allows for the detection and quantification of endogenous levels of natural and/or recombinant Human Eotaxin-2 proteins within the range of 8-500 pg/ml.

Manual (PDF): [Download](#)

NCBI Gene Symbol: CCL24

Specificity: The LumiAb<sup>™</sup>, Human Eotaxin-2 ELISA is capable of recognizing both recombinant and naturally produced Human Eotaxin-2 proteins. The antigens listed below were tested at 50 ng/ml and did not exhibit significant cross reactivity or interference. □ Human: Eotaxin, Eotaxin-3, MIP-1 (66 aa), MIP-1 (70 aa), MCP-3 □ Murine: Eotaxin, Eotaxin-2

**Sub Type:** None

**Synonyms:** C-C motif chemokine 24, CK-beta-6, Eosinophil chemotactic protein 2, Eotaxin-2, Myeloid progenitor inhibitory factor 2, Small-inducible cytokine A24, MPIF-2

**TargetName:** Human Eotaxin-2