

## Assay Biotechnology Company, Inc.

1276 Hammerwood Ave.

Sunnyvale, CA 94089

Phone: (408) 747-0185

Fax: (408)747-0145

www.assaybiotech.com

## GABRD Antibody

Catalog Number: C15856

### Description

#### Details

Product Name: GABRD Antibody

Catalog No: C15856

Stock:

#### Properties

Storage/Stability: Stable at -20°C for at least 1 year.

Applications: WB, ELISA

Clonality: Polyclonal

Host Species: Rabbit

Reactivity: Human, Mouse, Rat

Application Dilutions: WB: 1:500~1:1000 ELISA: 1:5000

Database Links: Gene ID: 2563, Swiss-Prot #: O14764, OMIM #: 137163/604233, Unigene #: Hs.113882

Immunogen: The antiserum was produced against synthesized peptide derived from human GABRD.

Immunogen Range: 141-190

MW: 50 kDa

NCBI Gene Symbol: GBRD

Physical Form: Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4,

150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Purification:** The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.

**Specificity:** GABRD Antibody detects endogenous levels of total GABRD protein.

**Synonyms:** Gamma-aminobutyric acid receptor subunit delta, GABA(A) receptor subunit delta

**Target Modification:** Unmodified/Total