

Assay Biotechnology Company, Inc.

1276 Hammerwood Ave.

Sunnyvale, CA 94089

Phone: (408) 747-0185

Fax: (408)747-0145

www.assaybiotech.com

Dynamin-1 Colorimetric Cell-Based ELISA Kit

Catalog Number: CB5212

Description

Details

Product Name: Dynamin-1 Colorimetric Cell-Based ELISA Kit

Catalog No: CB5212

Stock:

Properties

Detection Method: Colorimetric 450 nm

Dynamic Range: > 5000 Cells

Storage/Stability: 4°C/6 Months

Reactivity: Human, Mouse, Rat

Assay Type: Cell-Based

Database Links: Gene ID: 1759, UniProt ID: Q05193, OMIM #: 602377, Unigene #: Hs.522413

Format: 96-Well Microplate

Manual (PDF): [Download](#)

NCBI Gene
Symbol: DNM1

Sub Type: None

Synonyms: B-dynamin, D100, DNM, DNM1, Dynamin BREDNM19, Dynamin brain, dynamin-1

TargetName: Dynamin-1