

## Assay Biotechnology Company, Inc.

1276 Hammerwood Ave.

Sunnyvale, CA 94089

Phone: (408) 747-0185

Fax: (408)747-0145

www.assaybiotech.com

## EMR1 Antibody

Catalog Number: G090

## Description

### Details

Product Name: EMR1 Antibody

Catalog No: G090

Stock:

### Properties

Storage/Stability: Stable at -20°C for at least 1 year.

Applications: WB, ELISA

Clonality: Polyclonal

Host Species: Rabbit

Reactivity: Human

Application Dilutions: WB: 1:500~1:1000 ELISA: 1:1000

Database Links: Gene ID: 2015, Swiss-Prot #: Q14246, OMIM #: 600493, Unigene #: Hs.2375/

Immunogen: The antiserum was produced against synthesized peptide derived from human EMR1.

Immunogen Range: 818-867

MW: 97 kDa

NCBI Gene Symbol: EMR1

Physical Form: Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Purification:** The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.

**Specificity:** EMR1 Antibody detects endogenous levels of total EMR1 protein.

**Synonyms:** EGF-like module-containing mucin-like hormone receptor-like 1, Cell surface glycoprotein EMR1, EMR1 hormone receptor, EMR1, TM7LN3

**Target Modification:** Unmodified/Total