

Assay Biotechnology Company, Inc.

1276 Hammerwood Ave.

Sunnyvale, CA 94089

Phone: (408) 747-0185

Fax: (408)747-0145

www.assaybiotech.com

GPR108 Antibody

Catalog Number: G280

Description

Details

Product Name: GPR108 Antibody

Catalog No: G280

Stock:

Properties

Storage/Stability: Stable at -20°C for at least 1 year.

Applications: WB, IF, ELISA

Clonality: Polyclonal

Host Species: Rabbit

Reactivity: Human

Application Dilutions: WB: 1:500~1:1000 IF: 1:100~1:500 ELISA: 1:1000

Database Links: Gene ID: 56927, Swiss-Prot #: Q9NPR9, Unigene #: Hs.167641/

Immunogen: The antiserum was produced against synthesized peptide derived from human GPR108.

Immunogen Range: 155-204

MW: 60 kDa

NCBI Gene Symbol: GPR108

Physical Form: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Purification: The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.

Specificity: GPR108 Antibody detects endogenous levels of total GPR108 protein.

Synonyms: Protein GPR108, Lung seven transmembrane receptor 2, GPR108, LUSTR2

Target Modification: Unmodified/Total