

## Assay Biotechnology Company, Inc.

1276 Hammerwood Ave.

Sunnyvale, CA 94089

Phone: (408) 747-0185

Fax: (408)747-0145

www.assaybiotech.com

## GANP Colorimetric Cell-Based ELISA Kit

Catalog Number: CB5287

### Description

#### Details

Product Name: GANP Colorimetric Cell-Based ELISA Kit

Catalog No: CB5287

Stock:

#### Properties

Detection Method: Colorimetric 450 nm

Dynamic Range: > 5000 Cells

Storage/Stability: 4°C/6 Months

Reactivity: Human, Mouse

Assay Type: Cell-Based

Database Links: Gene ID: 8888, UniProt ID: O60318, OMIM #: 603294, Unigene #: Hs.389037

Format: 96-Well Microplate

Manual (PDF): [Download](#)

NCBI Gene  
Symbol: MCM3AP

Sub Type: None

Synonyms: 80 kda MCM3-associated protein, GANP protein, KIAA0572, MAP80, MCM3A, MCM3AP

TargetName: GANP